

Answers to UCLSE New Curriculum Chemistry



Experiment 6
A high school science fair showcased various experiments, including one investigating the rate of hydrogen gas production when magnesium reacts with two different solutions of dilute hydrochloric acid, labelled C and D, each with varying concentrations. Participants were recommended to use the set-up below.

Task:
(a) As a participant in the science fair,

Let C - more concentrated (1.0M)
D - less concentrated (0.7M)

(i) design an experiment to investigate the rate at which hydrogen gas is made when magnesium reacts with two different solutions of dilute hydrochloric acid, C and D, with different concentrations

Aim: To investigate the rate of evolution of Hydrogen gas when Magnesium ribbon reacts with two different acids of different concentration.

Variables: (i) Dependent: Volume of hydrogen gas produced.
(ii) Independent: Time taken during experiment
(iii) Controlled: Mass of magnesium used and volume of acid added.

Hypothesis: The rate of evolution of hydrogen gas depends on the high concentration of the acid solution. (More concentrated acid produces high rate of hydrogen gas).

Risks:

- Acid spilling on skin
- Cut from broken glass
- Choking fumes

Mitigation:

- Be careful and clean the table
- wear gloves
- use face masks

Practical Workbook

WAKATA



ANSWERS TO UCLSE NEW CURRICULUM CHEMISTRY PRACTICAL WORK BOOK



WAKATA PUBLISHERS

Answers to UCLSE New Curriculum Chemistry Practical Workbook

Sample copy

BY



WAKATA

TEL: +256 200 905 486 /+256 702 019 043 / +256 782 685 163

LOCATION: Bat Valley Complex-Bombo Rd. Kampala

WEBSITE: <http://wakatastationers.com>

EMAIL: wakatabureau@gmail.com

FACE BOOK: Wakiso Kampala Teachers' Association

WhatsApp: +256 782685163/+256 702 019 043

Experiment 1

A chemistry teacher conducted an experiment to investigate the rate of reaction between dilute hydrochloric acid and aqueous sodium thiosulphate. The experiment involved varying the volumes of the reactants to observe their impact on the rate of reaction.


Sodium thiosulphate reacts with hydrochloric acid according to the following equation.



Task:

(a) As a student conducting the experiment;

- (i) design an experiment to measure the rate of reaction between dilute hydrochloric acid and aqueous sodium thiosulphate. The setup below may be useful to you.



Bases of assessment:

Aim: To investigate /measure the reaction between $Na_2S_2O_3$ and HCl at different concentrations.

variables:

- (i) Dependent variable: Time taken for the cross to disappear.
(ii) Independent variable: Volume of sodium sulphate used.
(iii) Controlled variable: Volume of acid (BAZ) added.

Hypothesis: The rate of chemical reaction between $\text{Na}_2\text{S}_2\text{O}_3$ and HCl increases with increase in the concentration of the thiosulphate.

OR: The rate of chemical reaction between $\text{Na}_2\text{S}_2\text{O}_3$ and HCl decreases with increase in dilution or with decrease in the concentration

Apparatus:

- Stop clock or stop watch
- Beakers
- Measuring cylinder
- Conical flask

Risks:

- Spilling solution on table Be careful
- Acids can burn Use / wear protective gears

Procedure:

- Mark a cross (x) with a blue or black pen on the filter paper.
- Using a measuring cylinder, measure about 50cm^3 of BA1 (Sodium thiosulphate into a clean conical flask.
- Add 10cm^3 of HCl (BA2) to BA1 in the conical flask and at the same time start the stop clock.
- Shake/swirl the conical flask gently for the solution to mix and place it on the filter paper over the cross (x).
- Watch closely through the solution from above the conical flask. Stop the stop clock when the cross disappears and note the time taken for the cross to disappear.
- Pour away the mixture and thoroughly rinse the conical flask.
- Then measure 40cm^3 of BA1 followed by 10cm^3 of water into the conical flask.
- Add 10cm^3 of BA2 and at the same time start the stop clock.
- Follow the procedures from the first one again.
- Repeat the procedure with other 3 conical flasks and by measuring 30 , 20 and 10cm^3 of BA1 followed by 20 , 30 , 40cm^3 of water. The water should be added to BA1 before adding BA2. In each case keep the volume of BA2 the same.

(ii) carry out the experiment and record your findings.

Flask No.	1	2	3	4	5
Volume of BA1 (cm ³)	50	40	30	20	10
Volume of water (cm ³)	0	10	20	30	40
Volume of BA2 (cm ³)	10	10	10	10	10
Time, t in seconds	17.0	24.0	30.0	64.0	106.0
$\frac{1}{t} (s^{-1})$	0.059	0.042	0.033	0.0156	0.009

Sample copy

(iii) Describe the appearance of the mixture in the conical flask at the end of each experiment.

The mixture turns yellow. This is because the reaction between sodium thiosulphate and hydrochloric acid produces sulphur (colloidal) which covers the cross during the reaction.



- (iv) Explain, in terms of particles, why the rate of reaction was greatest in this experiment.

The rate of reaction is directly proportional to initial concentration of the reactant. The higher the concentration, the faster the rate of chemical reaction because increase in concentration increases the frequency of collision (the number of colliding / reacting particles in a given volume as a result the collisions between the particles increases).

- (b) Give the name of a more accurate piece of apparatus for measuring volumes than a measuring cylinder.

Burette, Pipette

- (c) Suggest the effect on the results of using a 100 cm^3 conical flask instead of a 250 cm^3 conical flask. Explain your answer.


No change on the results because the volume of a container, conical flask has no effect on the experimental result.

A graph of $\frac{1}{t}$ against volume of BA1 used



Experiment 6

A high school science fair showcased various experiments, including one investigating the rate of hydrogen gas production when magnesium reacts with two different solutions of dilute hydrochloric acid, labeled **C** and **D**, each with varying concentrations. Participants were recommended to use the set up below.



Task:

(a) As a participant in the science fair;

- (i) design an experiment to investigate the rate at which hydrogen gas is made when magnesium reacts with two different solutions of dilute hydrochloric acid, **C** and **D**, with different concentrations

Aim: To investigate the rate of evolution of Hydrogen gas when Magnesium ribbon reacts with two different acids of different concentration.

variables: **Dependent:** volume of hydrogen gas produced.
(ii) Independent: Time taken during experiment
(iii) Controlled: Mass of magnesium used and volume of acid added.

Hypothesis: The rate of evolution of hydrogen gas depends on the high concentration of the acid solution. (More concentrated acid produces high rate of hydrogen gas).

Risk:

- Acid spilling on skin
- Cut from broken glass
- Choking fumes

Mitigation

- Be careful and clean the table
- wear gloves
- use face masks

Procedure:

- A known mass of magnesium was cut and placed in a test tube / conical flask with a thread.
- 35cm³ of acid C was measured using a measuring cylinder and added to magnesium ribbon in a conical flask /test tube. At the same time a stop clock was started.
- The conical flask was connected to a gas syringe which is used to measure the volume of gas produced.
- The volume of hydrogen in the syringe is recorded at regular intervals until the reaction is completed. (When the whole mass of the magnesium ribbon has reacted).
- The above procedures were repeated with different acid D of different concentration (Less dilute).
- The results obtained are used to plot graphs of volume of hydrogen gas evolved against time on the same axes.

Sample copy

(ii) Conduct the experiment and record your findings

Length of Ribbon (cm)	1.0	2.0	3.0	4.0	5.0
Volume of $H_2(g)$ (cm^3) for solution C	3.0	5.4	7.5	9.8	10.8
Volume of $H_2(g)$ (cm^3) for solution D	2.2	3.6	5.2	9.2	10.8

OR:

Time (s)	0	30	60	120	150	180
Volume of $H_2(g)$ (cm^3) for solution C	0	3.0	4	9.8	10.8	10.8
Volume of $H_2(g)$ (cm^3) for solution D	0	2.2	3.6	5.2	9.2	10.8

(b) Explain what can be deduced about the concentrations of dilute hydrochloric acid C and dilute hydrochloric acid D.


The concentration of C is higher than that of D, the rate of reaction (production) of hydrogen gas is more with C than with D.

(c) Suggest why it is important to replace the bung in the conical flask immediately after adding the magnesium ribbon.

To avoid the gas from escaping

It is done to avoid wastage of gas produced.

A graph of volume of gas produced against time.



Experiment 35

A local municipality organized a workshop to educate residents on household cleaning products. During the workshop, participants conducted an experiment to observe the heat changes during the neutralization of acidic household waste with a basic solution.

As participants added a solution of acidic household waste to a solution of sodium hydroxide, they observed a temperature increase in the mixture. One participant noticed the temperature change and questioned why the household waste became warmer during the neutralization process.

The acidic solution provided is labeled **BA1** and the base provided is labeled **BA2**.

Task:

(a) As a learner of chemistry:

- Design an experiment to measure the temperature change and determine the amount of heat produced during the neutralization of **BA1** with **BA2**.

Aim: To measure the amount of heat and heat change produced during the neutralization of acidic household waste and sodium hydroxide.

Variables:

(i) Dependent:	Temperature change / rise
(ii) Independent:	Volume of acid added to NaOH
(iii) Controlled:	Volume of NaOH solution measured

Hypothesis: The reaction between **BA1** and **BA2** produces heat.

Risk:

Mitigation

- | | |
|----------------------------|--|
| - Acid spilling on skin | - Be careful and clean the working place |
| - Bad fumes from the waste | - wear face masks and other protective gears |

Procedures of experiment:

- 30cm^3 of BA2 is measured into a plastic beaker and its initial temperature noted and recorded as $T_1^\circ\text{C}$.
- The thermometer is rinsed and dried, then dipped into BA1 to obtain its initial temperature $T_2^\circ\text{C}$.
- Using a measuring cylinder, 30cm^3 of BA1 is measured and added to BA2 in plastic beaker.
- The solution is stirred with a thermometer and its maximum temperature rise is obtained, noted and recorded as $T_3^\circ\text{C}$.
- All the results are recorded.

Sample copy

- (ii) Conduct the experiment by mixing BA1 with BA2 and recording the temperature change over time.

- Initial temperature of BA1 = 25.0°C

- Initial temperature of BA2 = 26.0°C

- Average initial temperature = $\frac{25+26}{2}$

$$\equiv 25.5^{\circ}\text{C}$$

Volume of BA1 used = 30 cm³

Volume BA2 used = 30 cm³

Total volume of mixture = (30 + 30) cm³
= 60 cm³

Highest temperature attained by the mixture = 32.0 cm³

Sample copy

(iii) Calculate the maximum heat produced during the neutralization reaction.

$$\text{Temperature change } \theta = (\text{Final} - \text{Initial}) \text{ temperature}$$

$$= (32.0 - 25.5) {}^\circ\text{C}$$

$$\theta = 6.5 {}^\circ\text{C}$$

$$\text{Mass of the solution} = \text{Density} \times \text{volume}$$

$$1\text{g/cm}^3 \times 60\text{cm}^3$$

$$\text{Heat produced} = mc\theta$$

$$= 60 \times 4.2 \times 6.5$$

$$= 16385 \text{joules}$$

$$\text{Heat produced during the reaction} = -16385 \text{J}$$

(b) Based on your findings, explain to the participant how the neutralization reaction between BA₁ and BA₂ generates heat and helps in waste management.

- When BA₁ (acidic house hold waste) reacts with BA₂ (Sodium hydroxide) it produces -16385 J of heat implying the reaction is exothermic.
- Therefore during the disposal of waste products, neutralization has to be done to remove the acidic content in the waste before discharged.
- used in water treatment in order to reduce the damage of effluent.

Experiment Requirements (Practical Instructions)

Each candidate will require the following materials and apparatus.

Experiment 1

- (a) 250 cm³ of aqueous sodium thiosulfate, Na₂S₂O₃.5H₂O, containing 40 g / dm³ labelled **aqueous sodium thiosulfate**.
The aqueous sodium thiosulphate must be freshly prepared.
- (b) 100 cm³ of hydrochloric acid of concentration 2.0 mol / dm³ labelled **dilute hydrochloric acid**
- (c) access to water and distilled water
- (d) 50 cm³ measuring cylinder
- (e) 10 cm³ measuring cylinder
- (f) 250 cm³ conical flask
- (g) teat pipette
- (h) stop-clock or timer which can measure to an accuracy of 1 s

Per five candidates

A bucket labelled **quenching bath** must be provided.
The bucket must contain 1 dm³ of approximately 5% sodium carbonate solution (made up by dissolving 50 g of Na₂CO₃ or 135 g of Na₂CO₃.10H₂O in 1 dm³ of water) and Universal Indicator.
The Supervisor must monitor the colour of the Universal Indicator in each quenching bath to check that the solution has not become acidic. If the solution becomes acidic, the Supervisor must add more 5% sodium carbonate solution to the quenching bath.

Experiment 2

- Measuring cylinders (50cm³ and 10 cm³)
- conical flask
- filter paper
- stop clock (stop watch).

Experiment 3

- 50 cm³ of aqueous barium nitrate, Ba(NO₃)₂, of concentration 0.33 mol / dm³ labelled **aqueous barium nitrate**
- 50 cm³ of aqueous sodium carbonate, Na₂CO₃, of concentration 0.33 mol / dm³ labelled **aqueous sodium carbonate**
- access to water and distilled water
- 10 cm³ measuring cylinder

- 50 cm³ burette with stand and clamp
- funnel for filling burette
- 6 × **identical** test-tubes, capable of holding 15 cm³
- test-tube rack
- glass stirring rod
- stop-clock or timer which can measure to an accuracy of 1 s
- marker pen to write on glass
- ruler with millimetre graduations
- teat pipe

Experiment 4

- 50 cm³ of Sodium hydroxide solution (BA2) - 2M
- 50 cm³ of Sulphuric acid solution (BA1) - 1M
- Burette
- Pipette
- Conical flask
- Phenolphthalein indicator
- Distilled water
- Funnel
- Clamp stand
- White tile or paper

Experiment 5

- 200 cm³ of hydrochloric acid of concentration 0.20mol / dm³ in a stoppered container
- 100 cm³ of **solution E**, which is aqueous sodium carbonate of concentration 0.10 mol / dm³ in a stoppered container
- 50 cm³ of **solution F**, which aqueous sodium carbonate of concentration 0.15 mol / dm³ in a stoppered container
- thymolphthalein indicator
- methyl orange indicator
- 50 cm³ burette with clamp and stand
- white tile
- funnel to fill burette
- 25 cm³ measuring cylinder
- 250 cm³ conical flask
- dropping pipettes



- access to water and distilled water

Experiment 6

- 70 cm³ of dilute hydrochloric acid **C** of concentration 1.0 mol / dm³
- 70 cm³ of dilute hydrochloric **D** acid of concentration 0.7 mol / dm³
- Two 5 cm lengths of magnesium ribbon. This should be cleaned with abrasive paper and then coiled into a spiral by wrapping it around a pencil or similar to facilitate dropping it into a flask.
- 50 cm³ measuring cylinder
- stop-clock or timer which can measure to an accuracy of 1 s
- stirring thermometer, -10 °C to +110 °C with 1 °C graduations
- access to water and distilled water
- dropping pipette

Experiment 7

- 50 cm³ of solution A aqueous sodium hydroxide of concentration 0.10 mol / dm³
- 50 cm³ of solution B aqueous sodium hydroxide of concentration 0.05 mol / dm³
- 100 cm³ of dilute hydrochloric acid of concentration 0.15 mol / dm³
- thymolphthalein indicator
- access to water and distilled water
- 25 cm³ measuring cylinder
- 250 cm³ conical flask
- 150 cm³ burette stand and clamp for burette
- white tile
- small funnel to fill burette
- dropping pipettes

Experiment 8

- 4 g of anhydrous sodium carbonate, Na₂CO₃, in a dry stoppered test-tube, labelled **solid N**
The anhydrous sodium carbonate should be newly purchased.
- 6 g of hydrated sodium carbonate, Na₂CO₃•H₂O, in a stoppered test-tube, labelled **solid O**
Washing soda crystals are NOT suitable.
- access to water and at least 100 cm³ of distilled water
- 50 cm³ measuring cylinder

- 250 cm³ beaker for support
- 2 × polystyrene cups (cardboard cups may be used as an alternative)
- stop-clock or timer which can measure to an accuracy of 1 s
- stirring thermometer, -10 °C to 110 °C at 1 °C graduations
- teat pipette
- spatula

Experiment 9

- burette
- stand and clamp for burette
- white tile
- small funnel to fill burette
- 100 cm³ of **solution K**, which is aqueous sodium carbonate of concentration 0.20 mol / dm³
- 100 cm³ of **solution L**, which is aqueous sodium carbonate of concentration 0.10 mol / dm³
- 100 cm³ of **acid M**, which is dilute hydrochloric acid of concentration 0.40 mol / dm³
- 100 cm³ of **acid N**, which is dilute hydrochloric acid of concentration 0.20 mol / dm³
- 25 cm³ measuring cylinder
- 250 cm³ conical flask
- methyl orange indicator
- access to water and distilled water
- dropping pipettes

Experiment 10

- 150 cm³ of dilute hydrochloric acid of concentration 1.5 mol / dm³
- 30 cm³ of **solution H**, which is aqueous sodium hydroxide of concentration 1.5 mol / dm³.
- 30 cm³ **solution H**, which is aqueous sodium hydroxide of concentration 0.8 mol / dm³
- burette
- stand and clamp for burette
- small funnel to fill burette
- 50 cm³ measuring cylinder
- stirring thermometer, -10 °C to +110 °C with 1 °C graduations
- 100 cm³ beaker
- access to water and distilled water
- dropping pipette



Experiment 11

- (a) three 10 cm³ measuring cylinders
- (b) one 50 cm³ burette
- (c) two boiling tubes
- (d) 100 cm³ of aqueous potassium peroxodisupphate, K₂S₂O₈, of concentration 0.02 mol / dm³, (5.4 g /dm³), labelled solution L
- (e) 50cm³ of aqueous potassium iodide, KI, of concentration 0.5mol/dm³, (83g/dm³), labelled aqueous potassium iodide. This solution should be freshly prepared.
- (f) 50 cm³ of aqueous sodium thiosulfate, Na₂S₂O₃.5H₂O, of concentration 0.01 mol / dm³, (2.5 g /dm³), labelled aqueous sodium thiosulphate
- (g) 20cm³ of 0.2% starch indicator solution, labelled starch indicator
- (h) access to water and distilled water
- (i) stop clock or timer
- (j) a Bunsen burner and matches
- (k) holder for heating the boiling tube
- (l) one stirring thermometer, 0–110 °C

Note: Mixing 10 cm³ of solution L with 5 cm³ of aqueous potassium iodide, 3 cm³ of aqueous sodium thiosulfate and 1 cm³ of starch solution should produce a blue colour within five minutes.

Experiment 12

- (a) one 25 cm³ or 50 cm³ measuring cylinder
- (b) one 50 cm³ burette with stand and clamp
- (c) one 250 cm³ conical flask
- (d) 100 cm³ of aqueous hydrochloric acid of concentration 0.2 mol /dm³, labelled hydrochloric acid
- (e) 50 cm³ of an aqueous solution which is a mixture of sodium hydroxide, NaOH, of concentration 0.05mol /dm³, and sodium carbonate, Na₂CO₃, of concentration 0.05 mol /dm³ labelled solution R

Solution R can be made up by dissolving 2 g of sodium hydroxide and 5.3 g of anhydrous sodium carbonate in distilled water and making the solution up to 1 dm³

- (f) access to water and distilled water
- (g) phenolphthalein indicator
- (h) methyl orange indicator

- (i) one white tile
- (j) one test-tube
- (k) teat pipette
- (l) one funnel for filling burette

Experiment 13

- (a) 50 cm³ measuring cylinder
- (b) stirring thermometer, 0–110 °C
- (c) polystyrene cup and 250 cm³ beaker
- (d) 150 cm³ of aqueous sodium hydroxide, NaOH, of concentration 1.0 mol / dm³
- (e) 100 cm³ of aqueous hydrochloric acid, HCl, of concentration 2.0 mol / dm³, labelled **acid A**
- (f) 100 cm³ of aqueous sulfuric acid, H₂SO₄, of concentration 2.0 mol / dm³, labelled **acid B**
- (g) access to water and distilled water
- (h) 50 cm³ burette with stand and clamp
- (i) teat pipette
- (j) funnel for filling burette

Experiment 14

- (a) one 25cm³ measuring cylinder
- (b) one 50cm³ burette with stand and clamp
- (c) one 250cm³ conical flask
- (d) 100cm³ of an aqueous solution of sodium carbonate, Na₂CO₃, of concentration 0.05mol/dm³, labelled aqueous sodium carbonate (This can be made by adding 5.3g of anhydrous sodium carbonate to 500cm³ of distilled water. Dissolve the solid and then make the volume up to 1dm³ with distilled water.)
- (e) 100cm³ of hydrochloric acid of concentration 0.10mol/dm³, labelled **solution A**
- (f) 100cm³ of hydrochloric acid of concentration 0.20mol/dm³, labelled **solution B**
- (g) access to water and distilled water
- (h) thymolphthalein indicator
- (i) methyl orange indicator
- (j) a white tile
- (k) funnel for filling burette
- (l) teat pipette

Experiment 15

- (a) 2 × 25cm³ measuring cylinders
- (b) 1 × 10cm³ measuring cylinder
- (c) 1 × 50cm³ burette with stand and clamp
- (d) 1 × 250cm³ conical flask
- (e) 100cm³ of an aqueous solution of hydrated sodium thiosulfate, Na₂S₂O₃.5H₂O, of concentration 0.1mol/dm³ containing 25g per dm³, labelled aqueous sodium thiosulfate.
3cm³ of 1mol/dm³ aqueous sodium hydroxide should be added to the above solution to ensure that it is not acidic.
- (f) 50cm³ of an aqueous solution of potassium iodate, KIO₃, containing 6g per dm³, labelled **solution C**
- (g) 50cm³ of an aqueous solution of potassium iodate, KIO₃, containing 3g per dm³, labelled **solution D**
- (h) 50cm³ of sulfuric acid of concentration 1mol/dm³, labelled dilute sulfuric acid
- (i) access to water and distilled water
- (j) two samples in stoppered test-tubes of 1g solid potassium iodide, each labelled 1g of potassium iodide
- (k) 10cm³ of starch solution
- (l) white tile
- (m) funnel for filling burette

Experiment 16

- 200cm³ of dilute hydrochloric acid of concentration 2.0mol/dm³ labelled as **dilute hydrochloric acid**
- a tube or weighing boat containing 1g of sodium hydrogencarbonate labelled as **1g of sodium hydrogencarbonate**
- a tube or weighing boat containing 2g of sodium hydrogencarbonate labelled as **2g of sodium hydrogencarbonate**
- a tube or weighing boat containing 3g of sodium hydrogencarbonate labelled as **3g of sodium hydrogencarbonate**
- a tube or weighing boat containing 5g of sodium hydrogencarbonate labelled as **5g of sodium hydrogencarbonate**
- a tube or weighing boat containing 6g of sodium hydrogencarbonate labelled as **6g of sodium hydrogencarbonate**
- a tube or weighing boat containing 7g of sodium hydrogencarbonate labelled as **7g of sodium hydrogencarbonate**

- 250cm³ conical flask
- 25cm³ measuring cylinder
- stop-clock or timer which can measure to an accuracy of 1s 1 stirring thermometer, -10°C to +110°C at 1°C graduations
- dropping pipettes
- access to water and distilled water

Experiment 17

- (a) 3g of anhydrous sodium carbonate in a stoppered test-tube labelled **A**;
- (b) 4g of potassium hydrogen carbonate in a stoppered test-tube labelled **B**;
- (c) two plastic drinking cups;
- (d) 200cm³ of hydrochloric acid of concentration 2mol/dm³;
- (e) stirring thermometer, covering the range 0°C to 20°C graduated in degrees and clamped securely in such a way as to allow a plastic cup of fluid to be swirled underneath;
- (f) a timer, e.g. a stop clock, stop watch or wall clock with seconds hand;
- (g) one measuring cylinder, 0cm³ to 50cm³.


When solid **A** is added to 30cm³ of hydrochloric acid there should be a measurable temperature rise.

When solid **B** is added to 30cm³ of hydrochloric acid there should be a measurable temperature decrease

Experiment 18

- (a) about 100cm³ of aqueous hydrogen peroxide of concentration 5 volume (1.5% w/v), labelled solution **A**.
- (b) about 5g of manganese (IV) oxide powder,
- (c) distilled water,
- (d)
 - (i) two measuring cylinders, one to measure 20cm³ of liquid and one to collect 50cm³ of gas
 - (ii) one conical flask, 100cm³, with bung; rubber connector and delivery tube;
 - (iii) two clamps, one to support the conical flask and one to support the inverted measuring cylinder;
 - (iv) trough, or other suitable container, of water,

The students will need to set this apparatus as below.



- (e) a stop clock or access to a timer with a second hand;
- (f) one spatula;
1 spatula measure of manganese (IV) oxide added to 20cm³ of 5 volume hydrogen peroxide should produce atleast 20cm³ and not more than 40cm³ of gas in 20 seconds;
- (g) splints;
- (h) Bunsen burner;
- (i) test tube.

Experiment 19

- (a) one 50cm³ measuring cylinder;
- (b) one 50cm³ burette;
- (c) one or more 250cm³ conical flasks;
- (d) 50cm³ of aqueous sulphuric acid of concentration 2.0 mol/dm³, labelled acid C;
- (e) 50cm³ of aqueous nitric acid of concentration 2.0 mol/dm³, labelled acid D;
- (f) 50cm³ of aqueous sodium hydroxide of concentration 2.0mol/dm³, labelled aqueous sodium hydroxide;
- (g) access to water and distilled water;
- (h) one stirring thermometer, 0°C–110°C;
- (i) a white tile;
- (j) phenolphthalein indicator, labelled phenolphthalein indicator .

20cm³ of sodium hydroxide solution + 10 cm³ of acid C should give a temperature rise of about 10°C.

Experiment 20

- (a) Bunsen burner and matches;
- (b) 12g of potassium nitrate crystals in a stoppered boiling tube, labelled salt A;
- (c) one 50cm³ burette;
- (d) one 250cm³ beaker,
- (e) distilled water;
- (f) one stirring thermometer, covering the range 0°C to 100°C.

12g of potassium nitrate dissolved in 10cm³ of hot distilled water should form crystals on cooling at about 70°C.

Experiment 21

- Potassium chloride
- (a) Bunsen burner and matches;
- (b) 12g of potassium chloride crystals in a stoppered boiling tube, labelled potassium chloride crystals;
- (c) one 50cm³ burette;
- (d) one 250cm³ beaker,
- (e) distilled water;
- (f) one stirring thermometer, covering the range 0°C–100°C.

Experiment 22

- (a) 50g of Potassium chloride (KCl)
- (b) Distilled water
- (c) Thermometer
- (d) Stirring rod
- (e) Weighing balance
- (f) Beakers
- (g) Graduated cylinder
- (h) Hot plate or water bath
- (i) Stopwatch or timer

Experiment 23

- 50cm³ of BA1
- 50cm³ of BA2
- Thermometer 100cm³ measuring cylinder
- 250cm³ conical flask / titration flask
- 50cm³ burette
- 20.cm³ or 25.0cm³ pipette
- A stand
- 250ml plastic beakers
- BA2 is a 1M sodium hydroxide
- BA1 is a 1M sulphuric acid

Experiment 24

- 50cm³ of **BA1**
- 50cm³ of **K**
- Thermometer
- 100cm³ measuring cylinder
- 250cm³ conical flask / titration flask
- 50cm³ burette
- 20.cm³ or 25.0cm³ pipette
- A stand
- 2 250ml plastic beaker
- **BA1** is a 0.1M sodium hydroxide
- **K** is a 1M hydrochloric acid

Experiment 25

- 50cm³ of **BA1**
- 50cm³ of **BA2**
- Thermometer
- 100cm³ measuring cylinder
- 250cm³ conical flask / titration flask
- 50cm³ burette
- 20.cm³ or 25.0cm³ pipette
- A stand
- 2 250ml plastic beakers
- **BA2** is a 0.1M sodium carbonate
- **BA1** is a 1M sulphuric acid

Experiment 26

- 50cm³ of **BA1**
- 50cm³ of **BA2**
- 100cm³ measuring cylinder
- 250cm³ conical flask / titration flask
- 50cm³ burette
- 20.cm³ or 25.0cm³ pipette
- A stand
- 2 250ml plastic beakers
- Phenolphthalein and methyl orange indicators
- **BA2** is a 0.1M sodium hydroxide
- **BA1** is a 0.5M Hydrochloric acid

Experiment 27

- 50cm³ of **BA1**
- 50cm³ of **BA2**
- 100cm³ measuring cylinder
- 250cm³ conical flask / titration flask
- 50cm³ burette
- 20.cm³ or 25.0cm³ pipette
- A stand
- 2 250ml plastic beakers
- Phenolphthalein and methyl orange indicators
- **BA2** is a 0.1M sodium hydroxide
- **BA1** is a 1M sulphuric acid

Experiment 28

- 50cm³ of **BA1**
- 50cm³ of **BA2**
- 100cm³ measuring cylinder
- 250cm³ conical flask / titration flask
- 50cm³ burette
- 20.cm³ or 25.0cm³ pipette
- A stand
- 2 250ml plastic beakers
- Phenolphthalein and methyl orange indicators
- **BA2** is a 0.1M sodium hydroxide
- **BA1** is a 0.2M acetic acid

Experiment 29

- 50cm³ of **BA1**
- 50cm³ of **BA2**
- 100cm³ measuring cylinder
- 250cm³ conical flask / titration flask
- 50cm³ burette
- 20.cm³ or 25.0cm³ pipette
- A stand
- 2 250ml plastic beakers
- Phenolphthalein and methyl orange indicators
- **BA2** is a 0.1M sodium hydroxide
- **BA1** is a 0.3M citric acid

Experiment 30

- 50cm³ of **HCl**
- 50cm³ of **R**
- 100cm³ measuring cylinder
- 250cm³ conical flask / titration flask
- 50cm³ burette
- 20.cm³ or 25.0cm³ pipette
- A stand
- 2 250ml plastic beakers
- Phenolphthalein and methyl orange indicators
- **HCl** is a 0.1M Hydrochloric acid
- **R** is river water

Experiment 31

- 50cm³ of **KOH**
- 50cm³ of **J**
- 100cm³ measuring cylinder
- 250cm³ conical flask / titration flask
- 50cm³ burette
- 20.cm³ or 25.0cm³ pipette
- A stand
- 2 250ml plastic beakers
- Phenolphthalein and methyl orange indicators
- **KOH** is a 0.1M Potassium hydroxide
- **J** is orange Juice sample

Experiment 32

- 50cm³ of **BA1**
- 50cm³ of **BA2**
- 50cm³ measuring cylinder
- 250cm³ conical flask / titration flask
- 50cm³ burette
- 20.cm³ or 25.0cm³ pipette
- A stand
- 2 250ml plastic beakers
- Thermometer
- Phenolphthalein and methyl orange indicators
- **BA2** is a 0.1M sodium hydroxide
- **BA1** is a 0.1M Hydrochloric acid

Experiment 33

- 50cm³ of **BA1**
- 50cm³ of **BA2**
- 50cm³ measuring cylinder
- 250cm³ conical flask / titration flask
- 50cm³ burette
- 20.cm³ or 25.0cm³ pipette
- A stand
- 2 250ml plastic beakers
- Thermometer
- Calorimeter
- Phenolphthalein and methyl orange indicators
- **BA2** is a 0.1M sodium hydroxide
- **BA1** is Acidic water solution (Contaminated water sample)

Experiment 34

- 50cm³ of **BA1**
- 50cm³ of **BA2**
- 50cm³ measuring cylinder
- 250cm³ conical flask / titration flask
- 50cm³ burette
- 20.cm³ or 25.0cm³ pipette
- A stand
- 2 250ml plastic beakers
- Thermometer
- Calorimeter
- Phenolphthalein and methyl orange indicators
- **BA2** is a 0.1M sodium hydroxide
- **BA1** is fatty acid solution (Mixture of fatty acids from vegetable oils)

Experiment 35

- 50cm³ of **BA1**
- 50cm³ of **BA2**
- 50cm³ measuring cylinder
- 250cm³ conical flask / titration flask
- 50cm³ burette
- 20.cm³ or 25.0cm³ pipette
- A stand
- 2 250ml plastic beakers
- Thermometer
- Calorimeter
- Phenolphthalein and methyl orange indicators
- **BA2** is a 0.1M sodium hydroxide

- **BA1** is Acidic household waste solution
(Diluted household acidic waste)

Experiment 36

- 50cm³ of **BA1**
- 50cm³ of **BA2**
- 50cm³ measuring cylinder
- 250cm³ conical flask / titration flask
- 50cm³ burette
- 20.cm³ or 25.0cm³ pipette
- A stand
- 2 250ml plastic beakers
- Thermometer
- Calorimeter
- Phenolphthalein and methyl orange indicators
- **BA2** is a 0.1M sodium hydroxide
- **BA1** is Acidic skincare ingredients solution
(Mixture of acidic components from skincare products)

Experiment 37

- 50cm³ of **BA1**
- 200cm³ of **BA2**
- 50cm³ measuring cylinder
- 250cm³ conical flask / titration flask
- 50cm³ burette
- 20.cm³ pipette
- A stand
- 2 250ml plastic beakers
- Thermometer
- Phenolphthalein and methyl orange indicators
- **BA2** is a 0.1M sodium hydroxide
- **BA1** is 0.1M hydrochloric acid

END